



Recombinant Human CCN family member 2 (CCN2)

Product Code	CSB-YP006147HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P29279
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	QNCS GPCRCPDEPA PRCPAGVSLV LDGCGCCRVC AKQLGELCTE RDPCDPHKGL FCHFGSPANR KIGVCTAKDG APCIFGGTVY RSGESFQSSC KYQCTCLDGA VGC MPLCSMD VRLPSPDCPF PRRVKLPGKC CEEWCDEPK DQTVVG PALA AYRLEDTFGP DPTMIRANCL VQTTEWSACS KTCGMGISTR VTNDNASCR L EKQSRLCMVR PCEADLEENI KKGKKCIRTP KISKPIKFEL SGCTSMKTYR AKFCGVCTDG RCCTPHRTTT LPVEFKCPDG EVMKKNMMFI KTCACHYNCP GDNDIFESLY YRKMYGDMA
Source	Yeast
Target Names	CCN2
Protein Names	Recommended name: Connective tissue growth factor Alternative name(s): CCN family member 2 Hypertrophic chondrocyte-specific protein 24 Insulin-like growth factor-binding protein 8 Short name= IBP-8 Short name= IGF-binding pr
Expression Region	27-349
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.