



Recombinant *Saccharomyces cerevisiae* Eukaryotic translation initiation factor 4B (TIF3)

Product Code	CSB-MP007554SVG
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P34167
Product Type	Recombinant Protein
Immunogen Species	<i>Saccharomyces cerevisiae</i> (strain ATCC 204508 / S288c) (Baker's yeast)
Purity	>85% (SDS-PAGE)
Sequence	MAPPKKTVKK MDLTSFLNDD TFGSSWAEED VDLNKITIP I ETANANTIPL SELAHAKNNS NNTRSGGF GG SFGGRSRLDP ALGGGSSDRR EEYPVPDAPP YRAVINNIPW DITPEGVQAW VEDGLVKPEA VEEVVLPKNL RDPTRLKGNA FVTLKERADL VAVLKFNGTK LNERTVYVSV AAPRRGGGAD VDWSSARGSN FQGDGREDAP DLDWGAARGS NFRGPRRERE EVDIDWTAAR GSNFQGSSRP PRREREEVDI DWSAARGSNF QGSSRPPRRE REEPDIDWSA ARGSNFQSSS RPPRRERE EP DIDWSAARGS NFQSSSRPPR REREKEEPAL DWGAARGAQF GKPQQTKNTY KDRSLTNKKT TDEQPKIQKS VYDVLRTEDD DEDEEA EKQN GDAKENKVDA AVEKLQDKTA QLTVEDGDNW EVLGKK
Source	Mammalian cell
Target Names	TIF3
Protein Names	Recommended name: Eukaryotic translation initiation factor 4B Short name=eIF-4B
Expression Region	1-436
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.