



Recombinant *Rhizopus niveus* Phosphoglycerate kinase 1 (PGK1)

Product Code	CSB-BP017856RDC
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P29405
Product Type	Recombinant Protein
Immunogen Species	<i>Rhizopus niveus</i>
Purity	>85% (SDS-PAGE)
Sequence	MSLSNKLSIR DLNLKDKRVL IRVDFNVPMK DGAITNNNRI VQALPTVKYA LDNGASAVIL MSHLGRPNGE AVAKYSLKPV AAEVEKLLGK PVEFLNDCVG PDVEKACQSA KDGKVILLEN LRFHIEEEGS AKVDGQKVKA DAEAIKKFRA SLTTLADIYI NDAFGTAHRA HSSVMGV DLS QRAAGFLMQK ELEYFAKALE NPSRPFLAIL GGAKVSDKIQ LIENMLDKVN ALIICGGMAF TFKKTLDNVK IGKSLFDEPG SKLVQNLVKK AAEKNVKIVF PVDFITADKF APDASTGYAT DDDGI PDGWQ GLDCGERSNK LFREEILKSK TIVWNGPSGV FEFD AFSSGT KAVLDAVINA TKEGATTIII GGGDTATAAL KWGAEGQVSH ISTGGGASLE LLEGKELPGV TALSSKN
Source	Baculovirus
Target Names	PGK1
Protein Names	Recommended name: Phosphoglycerate kinase 1 EC= 2.7.2.3
Expression Region	1-417
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.