



Recombinant *Macropus eugenii* Transthyretin (TTR)

Product Code	CSB-EP025270MPE
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P42204
Product Type	Recombinant Protein
Immunogen Species	<i>Macropus eugenii</i> (Tammar wallaby) (<i>Notamacropus eugenii</i>)
Purity	>85% (SDS-PAGE)
Sequence	AAVHHEGEHS KCPLMVKVLD AVRGRPAVNV DVKVFKKTEE QTWELFAAGK TNDNGEIHLEL TTDDKFGGGL YKVEFDTISY WKALGVSPFH EYADVFTAN DAGHRHYTIA ALLSPYSFST TAIVSNPTE
Source	E.coli
Target Names	TTR
Protein Names	Recommended name: Transthyretin Alternative name(s): Prealbumin
Expression Region	21-149
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.