



Recombinant *Saccharomyces cerevisiae* Phosphoserine phosphatase (SER2)

Product Code	CSB-MP018938SVG
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P42941
Product Type	Recombinant Protein
Immunogen Species	<i>Saccharomyces cerevisiae</i> (strain ATCC 204508 / S288c) (Baker's yeast)
Purity	>85% (SDS-PAGE)
Sequence	MSK FVITCIA HGENLPKETI DQIAKEITES SAKDVSINGT KKLSARATDI FIEVAGSIVQ KDLKNKLTNV IDSHNDVDVI VSDNEYRQA KKL FVFD MDS TLIYQEVIEL IAAYAGVEEQ VHEITERAMN NELDFKESLR ERVKLLQGLQ VDTLYDEIKQ KLEVTKGVPE LCKFLHKKNC KLAVLSGGFI QFAGFIKDQL GLDFCKANLL EVDTDGKLTG KTLGPIVDGQ CKSETLLQLC NDYNVPVEAS CMVGDGGNDL PAMATAGFGI AWNAKPKVQK AAPCKLNTKS MTDILYILGY TDDEIYNRQ
Source	Mammalian cell
Target Names	SER2
Protein Names	Recommended name: Phosphoserine phosphatase Short name= PSP Short name= PSPase EC= 3.1.3.3 Alternative name(s): O-phosphoserine phosphohydrolase
Expression Region	1-309
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.