



Recombinant Arabidopsis thaliana Mevalonate kinase (At5g27450)

Product Code	CSB-BP015247DOA
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P46086
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	>85% (SDS-PAGE)
Sequence	MEVKARAPGK IILAGEHAVV HGSTAVAAAI DLYTYVTLRF PLPSAENNDR LTLQLKDISL EFSWSLARIK EAIPYDSSTL CRSTPASCSE ETLKSIIVLV EEQNLPKEKM WLSSGISTFL WLYTRIIGFN PATVVINSEL PYGSGLGSSA ALCVLTAAL LASSISEKTR GNGWSSLDET NLELLNKWAF EGEKIIHGKP SGIDNTVSAY GNMIFCSGE ITRLQSNMPL RMLITNTRVG RNTKALVSGV SQRAVRHPDA MKSVFNAVDS ISKELAAIIQ SKDETSVTEK EERIKELMEM NQGLLLSMGV SHSSIEAVIL TTVKHKLVSK LTGAGGGGCV LLLLPTGTVV DKVVEELESS GFQCFTALIG GNGAQICY
Source	Baculovirus
Target Names	At5g27450
Protein Names	Recommended name: Mevalonate kinase Short name= MK EC= 2.7.1.36
Expression Region	1-378
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.