



Recombinant Tarsius bancanus Involucrin (IVL)

Product Code	CSB-MP011922TLI
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P24711
Product Type	Recombinant Protein
Immunogen Species	Cephalopachus bancanus (Western tarsier) (Tarsius bancanus)
Purity	≥85% (SDS-PAGE)
Sequence	MSQQQTLPVT LPPALSQELL KTVPPPANTQ QDQMKQPTPS PAPCQKGPSE LPVEKHPAPV KQVPEQECEP QQQDHQEPEL QLGRKQQEPQ EQEVHPGKQQ QKPQEQAHL GKKQEPQGQE VHLGKQQQKT QEQEVHLGKQ QQELQEQEVH LEKQLQEPQE VHLEKQLQEQ EVHLEKQLQE PEPELNLGKQ QQEPQEQAAY LGKQQQELPE PQDPELHLGK QQQEPQEVEV QLEKQQEAQE QELHLGKQQQ ESQEQLHLR KLQQVPQEPQ DQELHLGKQQ QELQEQEVHL GKQLQEPQEQ ELHLGRQQQE LQEEEVHLGM KEEQLLKHVE QQEGQLEQQE GQLKQPVCIP TPGVQDIQP AQPPKGEVLL PTEKQQKQEV QWPLKQE
Source	Mammalian cell
Target Names	IVL
Protein Names	Recommended name: Involucrin
Expression Region	1-387
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.