



Recombinant *Xenopus laevis* Peroxisome proliferator-activated receptor alpha (ppara)

Product Code	CSB-MP018421XBE
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P37232
Product Type	Recombinant Protein
Immunogen Species	<i>Xenopus laevis</i> (African clawed frog)
Purity	≥85% (SDS-PAGE)
Sequence	MSTIMVDTNS ELCILTPLDE DDLESPLSGE FLQDIVDIQD ITQTIGDDGS TPFGASEHQF FGNSPGSIGS VSTDLTDTLS PASSPASITF PAASGSAEDA ACKSLNLECR VCSDKASGFH YGVHACEGCK GFFRRTIRLK LVYDR CERMC KIQKKNRNKC QYCRFEKCLN VGMSHNAIRF GRMPRSEKAK LKAEVLMCDQ DVKDSQMADL LSLARLIYDA YLKNFNMNKV KARAILTGKA SNPPFVIHDM ETLCMAEKT L VAKLVANGIQ NKEAEVRIFH CCQCTSVETV TELTEFAKSI PGFTELDLND QVTLLKYGVY EAMFAMLASV MNKDGMLVAY GNGFITREFL KSLRKPIGDM MEPKFEFAMK FNALELDDSD LSLFVAALIC CGDRPGLVNI PSIEKMQESI VHV LKLHLQS NHPDDSFLFP KLLQKMADLR QLVTEHAQLV QTIKKTETDA ALHPLLQEIY RDMY
Source	Mammalian cell
Target Names	ppara
Protein Names	Recommended name: Peroxisome proliferator-activated receptor alpha Short name= PPAR-alpha Alternative name(s): Nuclear receptor subfamily 1 group C member 1
Expression Region	1-474
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.