



Recombinant Thiobacillus ferrooxidans Citrate synthase (gltA)

Product Code	CSB-MP006031TOA
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P51045
Product Type	Recombinant Protein
Immunogen Species	Acidithiobacillus ferrooxidans (Thiobacillus ferrooxidans)
Purity	>85% (SDS-PAGE)
Sequence	MAEPNFAPGL EGVAATQSSI SNIDGAAGLL SYRGFAIADL AAHSSFEEVA LLLLDGVLPG AADLERFDHG LRAHRQVKYN VREIMKFMPV TGHPMDMLHC AVASLGMFYP QQELSDAERG NTLHLDAMAM RIIARMPATIV AMWEQMRFGN DPISRPDLS HAANFLYMLS GREPDAHTK ILDSCLILHA EHTINASTFS VLVTGSTLTN PYHVIGGAIG TLAGPLHGGA NQKVVEEMLEE ISSVQQVGAY LDRKMANKEK IWGFGHRIYK TRDPRAVILK GMMEDMASHG NLRHSSLFEI AIEVERQATE RLGAQGIHAN VDFYSGVLYH EMGIKADLFT PIFAMARSAG WLAHWREQLA DNRIFRPTQV YTGEQDRRYV PVAQRT
Source	Mammalian cell
Target Names	gltA
Protein Names	Recommended name: Citrate synthase EC= 2.3.3.1
Expression Region	1-386
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.