



# Recombinant Carassius auratus Glial fibrillary acidic protein (gfap)

<b>Product Code</b>	CSB-EP009369DQG-B
<b>Storage</b>	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
<b>Uniprot No.</b>	P48677
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Carassius auratus (Goldfish)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	<p>REVDRVMGLN DRFASYIEKV RFLEQQNKML VAELNQLRGK EPSRLGDIYQ  EELRELRQV DGLNAGKARL EIERDNLASD LATLKQRLQE ENALRQEAEN  NLNTFRQDVD EAALNRVQLE RKIDALQDEI SFLRKVHEEE MRQLQEQLIA  QQVHVDLDVS KPDLTALKE IRAQFEAMAT SNMQETEEWY RSKFADLTDA  AGRNAEALRQ AKQEANEYRR QIQGLTCDLE SLRGSNESLE RQLREMEERF  AIETAGYQDT VARLEDEIQM LKEEMARHLQ EYQDLLNVKL ALDIEIATYR  KLEGEESRI TVPVQNFTNL QFRDTS�DTK LTPEAHVKRS IIVRTVETRD  GEIIEKSTTE RKDLP</p>
<b>Source</b>	E.coli
<b>Target Names</b>	gfap
<b>Protein Names</b>	Recommended name: Glial fibrillary acidic protein Short name= GFAP
<b>Expression Region</b>	1-365
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.