


🕜 Tel: +1-301-363-4651 🛛 🖂 Email: cusabio@cusabio.com 🧿 Website: www.cusabio.com 🍙

VASP Antibody

Product Code	CSB-PA004424
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P50552
Immunogen	Synthesized peptide derived from Human VASP around the non- phosphorylation site of S238.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat,Monkey
Tested Applications	WB, IHC, IF, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, IF:1:200-1:1000, ELISA:1:20000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity- chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	VASP; Vasodilator-stimulated phosphoprotein; VASP
Immunogen Species	Homo sapiens (Human)
Target Names	VASP
Image	HeLa Western Blot analysis of Hela cells



1