



# Recombinant Human Scavenger receptor class B member 1 (SCARB1), partial, Baculovirus

<b>Product Code</b>	CSB-BP845139HU
<b>Abbreviation</b>	SCARB1
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q8WTV0
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	>85% (SDS-PAGE)(SDS-PAGE)
<b>Source</b>	Baculovirus
<b>Target Names</b>	SCARB1
<b>Protein Names</b>	Recommended name: Scavenger receptor class B member 1 Short name= SRB1 Alternative name(s): CD36 and LIMPII analogous 1 Short name= CLA-1 CD36 antigen-like 1 Collagen type I receptor, thrombospondin receptor-like 1 SR-BI CD_an
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	N-terminal His-tagged and C-terminal Myc-tagged N-terminal His-tagged Tag-Free
<b>Protein Length</b>	Partial
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.