



Recombinant Human Scavenger receptor class B member 1 (SCARB1), partial, Mammalian cell

| | |
|--------------------------|---|
| Product Code | CSB-MP845139HU |
| Abbreviation | SCARB1 |
| Storage | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |
| Uniprot No. | Q8WTV0 |
| Product Type | Recombinant Protein |
| Immunogen Species | Homo sapiens (Human) |
| Purity | ≥85% (SDS-PAGE) |
| Source | Mammalian cell |
| Target Names | SCARB1 |
| Protein Names | Recommended name: Scavenger receptor class B member 1 Short name= SRB1 Alternative name(s): CD36 and LIMPII analogous 1 Short name= CLA-1 CD36 antigen-like 1 Collagen type I receptor, thrombospondin receptor-like 1 SR-BI CD_an |
| Notes | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week. |
| Tag Info | N-terminal His-tagged and C-terminal Myc-tagged N-terminal His-tagged Tag-Free |
| Protein Length | Partial |
| Shelf Life | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |