



# Recombinant *Marinobacter aquaeolei* Na (+)-translocating NADH-quinone reductase subunit E, partial

<b>Product Code</b>	CSB-EP378875MNI1
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	A1U1Y9
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	<i>Marinobacter hydrocarbonoclasticus</i> (strain ATCC 700491 / DSM 11845 / VT8)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Source</b>	E.coli
<b>Target Names</b>	nqrE
<b>Protein Names</b>	Recommended name: Na(+)-translocating NADH-quinone reductase subunit E Short name= Na(+)-NQR subunit E Short name= Na(+)-translocating NQR subunit E EC= 1.6.5.- Alternative name(s): NQR complex subunit E NQR-1 subunit E
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Partial
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.