



Recombinant Escherichia coli O157:H7 Glutathione synthetase (gshB)

Product Code	CSB-BP009969EOD
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P58578
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli O157:H7
Purity	>85% (SDS-PAGE)
Sequence	MIKLGIVMDP IANINIKKDS SFAMLLEAQR RGYELHYMEM GDLYLINGEA RAHTRTLNVK QNYEEWFSFV GEQDLPLADL DVILMRKDPP FDTEFIYATY ILERAEEKGT LIVNKPQSLR DCNEKLFATAW FSDLTPETLV TRNKAQLKAF WEKHSIILK PLDGMGGASI FRVKEGDPNL GVIAETLTHE GTRYCMAQNY LPAIKDGDKR VLVVDGEPVP YCLARIPQGG ETRGNLAAGG RGEPRPLTES DWKIARQIGP TLKEKGLIFV GLDIIGDRLT EINVTSPTEI REIEAEFPVS ITGMLMDAIE ARLQQ
Source	Baculovirus
Target Names	gshB
Protein Names	Recommended name: Glutathione synthetase EC= 6.3.2.3 Alternative name(s): GSH synthetase Short name= GSH-S Short name= GSHase Glutathione synthase
Expression Region	1-315
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.