



Recombinant Geobacter sp. Glycerol-3-phosphate acyltransferase (plsY), partial

| | |
|--------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Code | CSB-EP498625GFP1-B |
| Storage | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |
| Uniprot No. | B9M1A8 |
| Product Type | Recombinant Protein |
| Immunogen Species | Geobacter daltonii (strain DSM 22248 / JCM 15807 / FRC-32) |
| Purity | >85% (SDS-PAGE) |
| Source | E.coli |
| Target Names | plsY |
| Protein Names | Recommended name: Glycerol-3-phosphate acyltransferase Alternative name(s): Acyl-PO4 G3P acyltransferase Acyl-phosphate--glycerol-3-phosphate acyltransferase G3P acyltransferase Short name= GPAT EC= 2.3.1.n3 Lysophosphati |
| Notes | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week. |
| Tag Info | Tag type will be determined during the manufacturing process. |
| Protein Length | Partial |
| Reconstitution | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference. |
| Shelf Life | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |