




Recombinant Human Neuropeptide FF receptor 2 (NPFFR2), partial,Yeast

Product Code	CSB-YP015983HU1
Abbreviation	NPFFR2
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	A0PJM9
Product Type	Recombinant Protein
Purity	≥ 90% as determined by SDS-PAGE.
Research Area	Neuroscience
Source	Yeast
Protein Names	Recommended name: Neuropeptide FF receptor 2 Alternative name(s): G-protein coupled receptor 74 G-protein coupled receptor HLWAR77 Neuropeptide G-protein coupled receptor
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal His-tagged Tag-Free
Mol. Weight	50.7kDa
Protein Length	Full Length

Image



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.