



Recombinant Colobus polykomos Cytochrome b-c1 complex subunit Rieske, mitochondrial (UQCRFS1), partial

Product Code	CSB-YP721306DVI1
Abbreviation	UQCRFS1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q69BK0
Product Type	Recombinant Protein
Immunogen Species	Colobus polykomos (Western black-and-white colobus monkey)
Purity	≥85% (SDS-PAGE)
Source	Yeast
Target Names	UQCRFS1
Protein Names	Recommended name: Cytochrome b-c1 complex subunit Rieske, mitochondrial EC= 1.10.2.2 Alternative name(s): Complex III subunit 5 Cytochrome b-c1 complex subunit 5 Rieske iron-sulfur protein Short name= RISP Ubiquinol-cytoc
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Partial
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.