



# Recombinant Episoriculus fumidus NADH-ubiquinone oxidoreductase chain 4L (MT-ND4L), partial

|                          |   |
|--------------------------|---|
| <b>Product Code</b>      | CSB-BP846463EDAR1   |
| <b>Abbreviation</b>      | MT-ND4L   |
| <b>Storage</b>           | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.   |
| <b>Uniprot No.</b>       | Q953I6  |
| <b>Product Type</b>      | Recombinant Protein   |
| <b>Immunogen Species</b> | Episoriculus fumidus (Taiwan brown-toothed shrew) (Soriculus fumidus)   |
| <b>Purity</b>            | >85% (SDS-PAGE)   |
| <b>Source</b>            | Baculovirus   |
| <b>Target Names</b>      | MT-ND4L   |
| <b>Protein Names</b>     | Recommended name: NADH-ubiquinone oxidoreductase chain 4L EC= 1.6.5.3<br>Alternative name(s): NADH dehydrogenase subunit 4L   |
| <b>Notes</b>             | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.   |
| <b>Tag Info</b>          | Tag type will be determined during the manufacturing process.   |
| <b>Protein Length</b>    | Partial   |
| <b>Reconstitution</b>    | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference. |
| <b>Shelf Life</b>        | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.   |