



Recombinant *Saccharomyces monacensis* Protein ORM1 (ORM1), partial

Product Code	CSB-BP856867SAH1
Abbreviation	ORM1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q96495
Product Type	Recombinant Protein
Immunogen Species	<i>Saccharomyces pastorianus</i> (strain ATCC 76670 / Carlsberg bottom yeast no.2 / CBS 1503 / CLIB 180 / NBRC 10610 / NRRL Y-1525) (Saaz-type lager yeast) (<i>Saccharomyces monacensis</i>)
Purity	>85% (SDS-PAGE)
Source	Baculovirus
Target Names	ORM1
Protein Names	Recommended name: Protein ORM1
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Partial
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.