



Recombinant *Pseudomonas aeruginosa* Diacylglycerol kinase (dgkA), partial

| | |
|--------------------------|---|
| Product Code | CSB-EP872707EZ1-B |
| Abbreviation | dgkA |
| Storage | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |
| Uniprot No. | Q9HY23 |
| Product Type | Recombinant Protein |
| Immunogen Species | <i>Pseudomonas aeruginosa</i> (strain ATCC 15692 / DSM 22644 / CIP 104116 / JCM 14847 / LMG 12228 / 1C / PRS 101 / PAO1) |
| Purity | >85% (SDS-PAGE) |
| Source | <i>E. coli</i> |
| Target Names | dgkA |
| Protein Names | Recommended name: Diacylglycerol kinase Short name= DAGK EC= 2.7.1.107 Alternative name(s): Diglyceride kinase Short name= DGK |
| Notes | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week. |
| Tag Info | Tag type will be determined during the manufacturing process. |
| Protein Length | Partial |
| Reconstitution | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference. |
| Shelf Life | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |