



Recombinant *Pyrococcus woesei* Glyceraldehyde-3-phosphate dehydrogenase (gap)

Product Code	CSB-MP009232FHZ
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P61880
Product Type	Recombinant Protein
Immunogen Species	<i>Pyrococcus woesei</i>
Purity	>85% (SDS-PAGE)
Sequence	MKIKVGINGY GTIGKRVAYA VTKQDDMELI GVTKTKPDFE AYRAKELGIP VYAASEEFLP RFEKAGFEVE GTLNDLLEKV DIIVDATPGG MGEKNKQLYE KAGVKAIFQG GEKAEVAQVS FVAQANYEAA LGKDYVRVVS CNTTGLVRTL NAIKDYVDYV YAVMIRRAAD PNDIKRGPIN AIKPSVTIPS HHGPDVQTVI PINIETSAFV VPTTIMHVHS IMVELKKPLT REDVIDIFEN TTRVLLFEKE KGFESTAQLI EFARDLHREW NNLYEIAVWK ESINVKGNRL FYIQAVHQES DVIPENIDAI RAMFEIAEKW ESIKKTNKSL GILK
Source	Mammalian cell
Target Names	gap
Protein Names	Recommended name: Glyceraldehyde-3-phosphate dehydrogenase Short name= GAPDH EC= 1.2.1.59 Alternative name(s): NAD(P)-dependent glyceraldehyde-3-phosphate dehydrogenase
Expression Region	1-334
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.