




Recombinant Human Cytochrome P450 4F8 (CYP4F8), partial

Product Code	CSB-EP006455HU1
Abbreviation	Recombinant Human CYP4F8 protein, partial
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P98187
Form	Liquid or Lyophilized powder
Storage Buffer	If the delivery form is liquid, the default storage buffer is Tris/PBS-based buffer, 5%-50% glycerol. If the delivery form is lyophilized powder, the buffer before lyophilization is Tris/PBS-based buffer, 6% Trehalose.
Product Type	Recombinant Human Cytochrome P450 4F8(CYP4F8) ,partial
Immunogen Species	Homo sapiens (Human)
Purity	≥ 85% as determined by SDS-PAGE.
Sequence	TYAFYHNGRRLRCFPQPRKQNWFLGHLGLVTPTEEGLRVLTQLVATYPQGFV RWLGPITPIINLCHPDIVRSVINTSDAITDKDIVFYKTLKPWLGDGLLSVGDKWR HHRLLTPAFHFNILKPYIKIFSKSANIMHAKWQRLAMEGSTCLDVFEHISLMTL DSLQKCIFSFDNSNCQEKPEYITAIMELSALVVKRNNQFFRYKDFLYFLTPCGR RFHRACRLVHDFDTDAVIQERRRLTSQGVDDFLQAKAKSKTLDLDFIDVLLSLEDK NGKELSDDIRAEADTFMFGGHDTTASGLSWVLYNLARHPEYQERCRQEVQE LLKDREPKEIEWDDLAQLPFLTMCLKESLRLHPPIPTFARGCTQDVVLPDSRVI PKGNCVNCINIFAIHHNPSVWPDPEVYDPFRFDPENAQKRSPMAFIPFSAGPRN CIGQKFAMAEMKVVLALTLLRFRILPDHREPRRTPEIVLRAEDGLWLRVEPLG
Research Area	Cell Biology
Source	E.coli
Target Names	CYP4F8
Expression Region	38-520aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal 10xHis-tagged and C-terminal Myc-tagged
Mol. Weight	63.4 kDa
Protein Length	Partial
Image	



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.