



Recombinant Human Peripherin-2 (PRPH2), partial

| | |
|--------------------------|--|
| Product Code | CSB-YP018775HU1 |
| Storage | Store at -20°C, for extended storage, conserve at -20°C or -80°C. |
| Uniprot No. | P23942 |
| Product Type | Recombinant Protein |
| Immunogen Species | Homo sapiens (Human) |
| Purity | ≥85% (SDS-PAGE) |
| Source | Yeast |
| Target Names | PRPH2 |
| Protein Names | Recommended name: Peripherin-2 Alternative name(s): Retinal degeneration slow protein Tetraspanin-22 Short name= Tspan-22 |
| Notes | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week. |
| Tag Info | His tagged |
| Protein Length | partial |
| Shelf Life | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |