



Recombinant Human Peripherin-2 (PRPH2), partial

Product Code	CSB-YP018775HU1
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P23942
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Source	Yeast
Target Names	PRPH2
Protein Names	Recommended name: Peripherin-2 Alternative name(s): Retinal degeneration slow protein Tetraspanin-22 Short name= Tspan-22
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	His tagged
Protein Length	partial
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.