



Recombinant Human 5,6-dihydroxyindole-2-carboxylic acid oxidase (TYRP1), partial, Yeast

Product Code	CSB-YP025399HU1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P17643
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Source	Yeast
Target Names	TYRP1
Protein Names	Recommended name: 5,6-dihydroxyindole-2-carboxylic acid oxidase Short name= DHICA oxidase EC= 1.14.18.-Alternative name(s): Catalase B Glycoprotein 75 Melanoma antigen gp75 Tyrosinase-related protein 1 Short name= TRP S
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal His-tagged Tag-Free
Protein Length	Partial
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.