



Recombinant Polynucleobacter sp. ATP synthase subunit b (atpF), partial

| | |
|--------------------------|---|
| Product Code | CSB-YP002358PYU1 |
| Storage | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |
| Uniprot No. | A4SUT0 |
| Product Type | Recombinant Protein |
| Immunogen Species | Polynucleobacter asymbioticus (strain DSM 18221 / CIP 109841 / QLW-P1DMWA-1) (Polynucleobacter necessarius subsp. asymbioticus) |
| Purity | >85% (SDS-PAGE) |
| Source | Yeast |
| Target Names | atpF |
| Protein Names | Recommended name: ATP synthase subunit b Alternative name(s): ATP synthase F(0) sector subunit b ATPase subunit I F-type ATPase subunit b Short name= F-ATPase subunit b |
| Notes | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week. |
| Tag Info | Tag type will be determined during the manufacturing process. |
| Protein Length | Partial |
| Reconstitution | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference. |
| Shelf Life | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |