



# Recombinant *Saccharomyces cerevisiae*

## Adenylate kinase 1 (ADK1)

|                          |                                                                                                                                                                                                                                                                                                                                                                                                             |
|--------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Product Code</b>      | CSB-YP001508SVG                                                                                                                                                                                                                                                                                                                                                                                             |
| <b>Storage</b>           | Store at -20°C, for extended storage, conserve at -20°C or -80°C.                                                                                                                                                                                                                                                                                                                                           |
| <b>Uniprot No.</b>       | P07170                                                                                                                                                                                                                                                                                                                                                                                                      |
| <b>Product Type</b>      | Recombinant Protein                                                                                                                                                                                                                                                                                                                                                                                         |
| <b>Immunogen Species</b> | <i>Saccharomyces cerevisiae</i> (strain ATCC 204508 / S288c) (Baker's yeast)                                                                                                                                                                                                                                                                                                                                |
| <b>Purity</b>            | ≥85% (SDS-PAGE)                                                                                                                                                                                                                                                                                                                                                                                             |
| <b>Sequence</b>          | SSESIRMV LIGPPGAGKG TQAPNLQERF HAAHLATGDM LRSQIAKGTQ<br>LGLEAKKIMD QGGLVSDDIM VNMKDELTN NPACKNGFIL DGFPRTIPQA<br>EKLDQMLKEQ GTPLEKAIEL KVDELLVAR ITGRLIHPAS GRSYHKIFNP<br>PKEDMKDDVT GEALVQRSDD NADALKKRLA AYHAQTEPIV DFYKKTGIWA<br>GVDASQPPAT VWADILNKLK KD                                                                                                                                                |
| <b>Source</b>            | Yeast                                                                                                                                                                                                                                                                                                                                                                                                       |
| <b>Target Names</b>      | ADK1                                                                                                                                                                                                                                                                                                                                                                                                        |
| <b>Protein Names</b>     | Recommended name: Adenylate kinase 1 Short name= AK 1 EC= 2.7.4.3<br>Alternative name(s): ATP-AMP transphosphorylase 1 Adenylate kinase cytosolic and mitochondrial                                                                                                                                                                                                                                         |
| <b>Expression Region</b> | 3-222                                                                                                                                                                                                                                                                                                                                                                                                       |
| <b>Notes</b>             | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.                                                                                                                                                                                                                                                                                                         |
| <b>Tag Info</b>          | Tag type will be determined during the manufacturing process.                                                                                                                                                                                                                                                                                                                                               |
| <b>Protein Length</b>    | Full Length of Mature Protein                                                                                                                                                                                                                                                                                                                                                                               |
| <b>Reconstitution</b>    | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference. |
| <b>Shelf Life</b>        | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.<br>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.                                                                                                                          |