



Recombinant *Saccharomyces cerevisiae* Eukaryotic translation initiation factor 4E (CDC33)

Product Code	CSB-BP007556SVG
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P07260
Product Type	Recombinant Protein
Immunogen Species	<i>Saccharomyces cerevisiae</i> (strain ATCC 204508 / S288c) (Baker's yeast)
Purity	≥85% (SDS-PAGE)
Sequence	MSVEEVSKKF EENVSVDTT ATPKTVLSDS AHFDVKHPLN TKWTLWYTKP AVDKSESWSL LLRPVTSFQT VEEFWAIQN IPEPHELPLK SDYHVFNRNDV RPEWEDEANA KGGKWSFQLR GKGADIDELW LRTLLAVIGE TIDEDDSQIN GVVLSIRKGG NKFALWTKSE DKEPLLRI GG KFKQVLKLT D DGHLEFFPHS SANGRHPQPS ITL
Source	Baculovirus
Target Names	CDC33
Protein Names	Recommended name: Eukaryotic translation initiation factor 4E Short name= eIF-4E Short name= eIF4E Alternative name(s): eIF-4F 25 kDa subunit mRNA cap-binding protein
Expression Region	1-213
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.