



CD99 Antibody

Product Code	CSB-PA004973GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P14209
Immunogen	Recombinant Human CD99 protein
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, WB, IHC
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide and 50% glycerol pH 7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	CD99 antigen, E2 antigen, Protein MIC2, T-cell surface glycoprotein E2
Immunogen Species	Homo sapiens (Human)
Target Names	CD99