



# CD99 Antibody

<b>Product Code</b>	CSB-PA004973GA01HU
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	P14209
<b>Immunogen</b>	Recombinant Human CD99 protein
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Human
<b>Tested Applications</b>	ELISA, WB, IHC
<b>Form</b>	Liquid
<b>Conjugate</b>	Non-conjugated
<b>Storage Buffer</b>	PBS with 0.02% sodium azide and 50% glycerol pH 7.3.
<b>Purification Method</b>	Antigen Affinity Purified
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Alias</b>	CD99 antigen, E2 antigen, Protein MIC2, T-cell surface glycoprotein E2
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Target Names</b>	CD99
<b>Usage</b>	For Research Use Only. Not for use in diagnostic or therapeutic procedures.