



Recombinant Emericella nidulans Dolichyl pyrophosphate Glc1Man9GlcNAc2 alpha-1,3-glucosyltransferase (alg8), partial

| | |
|--------------------------|---|
| Product Code | CSB-EP691444EKF-B |
| Abbreviation | alg8 |
| Storage | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |
| Uniprot No. | Q5AWM9 |
| Product Type | Recombinant Protein |
| Immunogen Species | Emericella nidulans (strain FGSC A4 / ATCC 38163 / CBS 112.46 / NRRL 194 / M139) (Aspergillus nidulans) |
| Purity | >85% (SDS-PAGE) |
| Source | E.coli |
| Target Names | alg8 |
| Protein Names | Recommended name: Dolichyl pyrophosphate Glc1Man9GlcNAc2 alpha-1,3-glucosyltransferase EC= 2.4.1.265 Alternative name(s): Asparagine-linked glycosylation protein 8 Dol-P-Glc:Glc(1)Man(9)GlcNAc(2)-PP-dolichyl alpha-1,3-glucosyltransferas |
| Notes | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week. |
| Tag Info | Tag type will be determined during the manufacturing process. |
| Protein Length | Partial |
| Reconstitution | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference. |
| Shelf Life | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |