



Recombinant Ashbya gossypii Dol-P-Glc:Glc (2)Man (9)GlcNAc (2)-PP-Dol alpha-1,2-glucosyltransferase (ALG10), partial

Product Code	CSB-YP741649DOT
Abbreviation	ALG10
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q75AQ8
Product Type	Recombinant Protein
Immunogen Species	Ashbya gossypii (strain ATCC 10895 / CBS 109.51 / FGSC 9923 / NRRL Y-1056) (Yeast) (Eremothecium gossypii)
Purity	>85% (SDS-PAGE)
Source	Yeast
Target Names	ALG10
Protein Names	Recommended name: Dol-P-Glc:Glc(2)Man(9)GlcNAc(2)-PP-Dol alpha-1,2-glucosyltransferase EC= 2.4.1.256 Alternative name(s): Alpha-1,2-glucosyltransferase ALG10-A Alpha-2-glucosyltransferase ALG10 Asparagine-linked glycosylation protei
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Partial
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.