



Recombinant Escherichia coli O6 Citrate synthase (gltA)

Product Code	CSB-BP006031EGX
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P0ABH8
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli O6:H1 (strain CFT073 / ATCC 700928 / UPEC)
Purity	>85% (SDS-PAGE)
Sequence	MADTKAKLTL NGDTAVELDV LKGTLGQDVI DIRTLGSKGV FTFDPGFTST ASCESKITFI DGDEGILLHR GFPIDQLATD SNYLEVCYIL LNGEKPTQEQ YDEFKTTVTR HTMIHEQITR LFHAFRRDSH PMAVMCGITG ALAAFYHDSL DVNNPRHREI AAFRLLSKMP TMAAMCYKYS IGQPFVYPRN DLSYAGNFLN MMFSTPCEPY EVNPILERAM DRILILHADH EQNASTSTVR TAGSSGANPF ACIAAGIASL WGPAHGGANE AALKMLEEIS SVKHIPEFVR RAKDKNDSFR LMGFGHRVYK NYDPRATVMR ETCHEVLKEL GTKDDLLEVA MELENIALND PYFIEKKLYP NVDFYSGIIL KAMGIPSSMF TVIFAMARTV GWIAHWSEM SDGMKIARPR QLYTGYEKRD FKSDIKR
Source	Baculovirus
Target Names	gltA
Protein Names	Recommended name: Citrate synthase EC= 2.3.3.1
Expression Region	1-427
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.