



Recombinant Drosophila melanogaster Adenosine 3'-phospho 5'-phosphosulfate transporter 1 (sll)

Product Code	CSB-EP896082DLU
Abbreviation	sll
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q9VEI3
Product Type	Recombinant Protein
Immunogen Species	Drosophila melanogaster (Fruit fly)
Purity	>85% (SDS-PAGE)
Sequence	<p>MYAYNKMGRVPELVICSFIVVTLLVIHFFSDLLRASLGGYYNQDVTLSQLVESQ NSDYAW FLKLLVNCFGYSCVFVPGFLIYKYVGRINYLERGNKTFLLHKAINMCITGNSGYD QLDAGT STADKDRPAASTAPKRTSSQEAVQLLWCFGGLMISYLTWGVLQEKIMTQNYL NFTGESAK FKDSQFLVFSNRLLAFLVALAYLQWQPSPVRHRAPLYKYSYASFNIMSAWFQ YEALKFV NFPTQVLAKSCKIIPVMLMGKIMSKAKYESYEYVTALLISLGMIFMSGSSDSSK ASGVT TLTGIFLLSMYMFDSFTANWQGSFLKSYGMTPLQMMCGVNLFSSIFTGASLS MQGGFMD SLAFATEHPKFVFDMMVLSVCSAVGQLFIYHTIDVFGPVVFTIIMTLRQAVAIMLS CFIY QHSISLLGIFGVLVFAIFLRVYCTQRLRAIRKRAEANKPKMAV</p>
Source	E.coli
Target Names	sll
Protein Names	Recommended name: Adenosine 3'-phospho 5'-phosphosulfate transporter 1 Alternative name(s): PAPS transporter 1 Protein slalom
Expression Region	1-465
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.