



Recombinant Mesocricetus auratus Glucagon (GCG)

| | |
|--------------------------|---|
| Product Code | CSB-MP009315MRG |
| Storage | Store at -20°C, for extended storage, conserve at -20°C or -80°C. |
| Uniprot No. | P01273 |
| Product Type | Recombinant Protein |
| Immunogen Species | Mesocricetus auratus (Golden hamster) |
| Purity | >85% (SDS-PAGE) |
| Sequence | HSLQDTEEKS RSFPASQTDP LEDPDQINED KRHSQGTFTS DYSKYLDSRR AQDFVQWLMN TKRNRNINIA |
| Source | Mammalian cell |
| Target Names | GCG |
| Protein Names | Recommended name: Glucagon Cleaved into the following 8 chains: 1. Glicentin 2. Glicentin-related polypeptide Short name= 3. GRPP 4. Oxyntomodulin Short name= 5. OXM Short name= 6. OXY 7. Glucagon 8. Glucagon-I |
| Expression Region | 21-89 |
| Notes | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week. |
| Tag Info | Tag type will be determined during the manufacturing process. |
| Protein Length | Cytoplasmic domain |
| Reconstitution | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference. |
| Shelf Life | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |