



Recombinant Human coronavirus HKU1 Replicase polyprotein 1a (1a), partial

Product Code	CSB-MP314074HIU
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P0C6U3
Product Type	Recombinant Protein
Immunogen Species	Human coronavirus HKU1 (isolate N1) (HCoV-HKU1)
Purity	≥85% (SDS-PAGE)
Source	Mammalian cell
Target Names	1a
Protein Names	Recommended name: Replicase polyprotein 1a Short name= pp1a Alternative name(s): ORF1a polyprotein Cleaved into the following 11 chains: 1. Non-structural protein 1 Short name= 2. nsp1 Alternative name(s): p28 Non-structu
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Partial
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.