



Recombinant Human coronavirus NL63 Replicase polyprotein 1ab (rep), partial

| | |
|--------------------------|---|
| Product Code | CSB-MP315168HIX |
| Storage | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |
| Uniprot No. | P0C6X5 |
| Product Type | Recombinant Protein |
| Immunogen Species | Human coronavirus NL63 (HCoV-NL63) |
| Purity | ≥85% (SDS-PAGE) |
| Source | Mammalian cell |
| Target Names | rep |
| Protein Names | Recommended name: Replicase polyprotein 1ab Short name= pp1ab Alternative name(s): ORF1ab polyprotein Cleaved into the following 15 chains: 1. Non-structural protein 1 Short name= 2. nsp1 Alternative name(s): p9 Non-struct |
| Notes | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week. |
| Tag Info | Tag type will be determined during the manufacturing process. |
| Protein Length | Partial |
| Reconstitution | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference. |
| Shelf Life | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |