



Recombinant Human Leukocyte cell-derived chemotaxin 1 (LECT1), partial, Yeast

Product Code	CSB-YP012854HU
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	O75829
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Source	Yeast
Target Names	CNMD
Protein Names	Recommended name: Leukocyte cell-derived chemotaxin 1 Cleaved into the following 2 chains: 1. Chondrosurfactant protein Short name= 2. CH-SP 3. Chondromodulin-1 Alternative name(s): Chondromodulin-I Short name= ChM-I
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal His-tagged Tag-Free
Protein Length	Full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.