



Recombinant Lymphocytic choriomeningitis virus Pre-glycoprotein polyprotein GP complex (GPC), partial, Yeast

Product Code	CSB-YP357830LKW
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P09991
Product Type	Recombinant Protein
Immunogen Species	Lymphocytic choriomeningitis virus (strain Armstrong) (LCMV)
Purity	≥85% (SDS-PAGE)
Source	Yeast
Target Names	GPC
Protein Names	Recommended name: Pre-glycoprotein polyprotein GP complex Cleaved into the following 3 chains: 1. Stable signal peptide Short name= 2. SSP 3. Glycoprotein G1 Short name= 4. GP1 5. Glycoprotein G2 Short name= 6. GP2
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal His-tagged Tag-Free
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.