



Recombinant Human Scavenger receptor class B member 1 (SCARB1), partial, Yeast

Product Code	CSB-YP845139HU
Abbreviation	SCARB1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8WTV0
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)(SDS-PAGE)
Source	Yeast
Target Names	SCARB1
Protein Names	Recommended name: Scavenger receptor class B member 1 Short name= SRB1 Alternative name(s): CD36 and LIMPII analogous 1 Short name= CLA-1 CD36 antigen-like 1 Collagen type I receptor, thrombospondin receptor-like 1 SR-BI CD_an
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal His-tagged Tag-Free
Protein Length	Extracellular domain
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.