



Recombinant Human Adenosine receptor A2a (ADORA2A), partial, Mammalian cell

Product Code	CSB-MP001376HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P29274
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	RIREFRQTFRKIIRSHVLRQQEPFKAAGTSARVLAAHGSDGEQVSLRLNGHPP GVWANGS APHPERRPNGYALGLVSGGSAQESQGNTGLPDVELLSHELKGVCPPEPGLDD PLAQDGAGVS
Source	Mammalian cell
Target Names	ADORA2A
Protein Names	Recommended name: Adenosine receptor A2a
Expression Region	291-412
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal His-tagged and C-terminal Myc-tagged N-terminal His-tagged Tag-Free
Protein Length	partial, Cytoplasmic domain
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.