



A1BG Antibody, HRP conjugated

Product Code	CSB-PA001001LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P04217
Immunogen	Recombinant Human Alpha-1B-glycoprotein protein (284-408AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Alpha-1B-glycoprotein (Alpha-1-B glycoprotein), A1BG
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	A1BG
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.