



ABCB7 Antibody, HRP conjugated

Product Code	CSB-PA001053LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O75027
Immunogen	Recombinant Human ATP-binding cassette sub-family B member 7, mitochondrial protein (23-140AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	ATP-binding cassette sub-family B member 7, mitochondrial (ATP-binding cassette transporter 7) (ABC transporter 7 protein), ABCB7, ABC7
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	ABCB7
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.