



# ADCY6 Antibody, HRP conjugated

<b>Product Code</b>	CSB-PA001342LB01HU
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	O43306
<b>Immunogen</b>	Recombinant Human Adenylate cyclase type 6 protein (9-118AA)
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Human
<b>Tested Applications</b>	ELISA
<b>Form</b>	Liquid
<b>Conjugate</b>	HRP
<b>Storage Buffer</b>	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
<b>Purification Method</b>	>95%, Protein G purified
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Alias</b>	Adenylate cyclase type 6 (EC 4.6.1.1) (ATP pyrophosphate-lyase 6) (Adenylate cyclase type VI) (Adenylyl cyclase 6) (Ca(2+)-inhibitible adenylyl cyclase), ADCY6, KIAA0422
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Research Area</b>	Signal Transduction
<b>Target Names</b>	ADCY6