



ADCY6 Antibody, FITC conjugated

| | |
|----------------------------|---|
| Product Code | CSB-PA001342LC01HU |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | O43306 |
| Immunogen | Recombinant Human Adenylate cyclase type 6 protein (9-118AA) |
| Raised In | Rabbit |
| Species Reactivity | Human |
| Form | Liquid |
| Conjugate | FITC |
| Storage Buffer | Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4 |
| Purification Method | >95%, Protein G purified |
| Isotype | IgG |
| Clonality | Polyclonal |
| Alias | Adenylate cyclase type 6 (EC 4.6.1.1) (ATP pyrophosphate-lyase 6) (Adenylate cyclase type VI) (Adenylyl cyclase 6) (Ca(2+)-inhibitible adenylyl cyclase), ADCY6, KIAA0422 |
| Immunogen Species | Homo sapiens (Human) |
| Research Area | Signal Transduction |
| Target Names | ADCY6 |
| Usage | For Research Use Only. Not for use in diagnostic or therapeutic procedures. |