



ADCY6 Antibody, Biotin conjugated

Product Code	CSB-PA001342LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O43306
Immunogen	Recombinant Human Adenylate cyclase type 6 protein (9-118AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Adenylate cyclase type 6 (EC 4.6.1.1) (ATP pyrophosphate-lyase 6) (Adenylate cyclase type VI) (Adenylyl cyclase 6) (Ca(2+)-inhibitable adenylyl cyclase), ADCY6, KIAA0422
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	ADCY6
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.