



GRK3 Antibody, HRP conjugated

Product Code	CSB-PA001395LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P35626
Immunogen	Recombinant Human Beta-adrenergic receptor kinase 2 protein (343-480AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Beta-adrenergic receptor kinase 2 (Beta-ARK-2) (EC 2.7.11.15) (G-protein-coupled receptor kinase 3), GRK3, ADRBK2 BARK2
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	GRK3