



GRK3 Antibody, Biotin conjugated

Product Code	CSB-PA001395LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P35626
Immunogen	Recombinant Human Beta-adrenergic receptor kinase 2 protein (343-480AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Beta-adrenergic receptor kinase 2 (Beta-ARK-2) (EC 2.7.11.15) (G-protein-coupled receptor kinase 3), GRK3, ADRBK2 BARK2
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	GRK3
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.