



# UBA1 Antibody, Biotin conjugated

Product Code	CSB-PA00149D0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P22314
Immunogen	Recombinant Human Ubiquitin-like modifier-activating enzyme 1 protein (730-881AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Ubiquitin-like modifier-activating enzyme 1 (EC 6.2.1.45) (Protein A1S9) (Ubiquitin-activating enzyme E1), UBA1, A1S9T UBE1
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	UBA1